

Title: Porcelain Lamina Restoration – two cases of the use

MUDr. Markéta Polednová

Private office, Kolín

The objective of this presentation is to demonstrate two indications of the use of the porcelain laminate veneers. Perfect face is very important for many patients but on the other hand they are anxiously afraid of the big intervention on their teeth. Sometimes we cannot achieve a satisfying effect neither by whitening nor making resin composite restoration. Porcelain laminate veneers restoration is minimally invasive method of the treatment with good results regardless of the shade of the prepared tooth.

Case 1

A 41-year-old nonsmoker woman, with an excellent dental hygiene had a problem with the middle left incisor which was rotated and there was the resin composite filling. Patient was not still satisfied with her smile so we resolved it by porcelain veneers on the both middle incisors which corrected profile of teeth. Both middle incisors were vital.



Case 2

A nonsmoker man, aged 37, with good dental hygiene had the right middle incisor for many years nonvital. There was good endodontic treatment but the shade of this tooth was much more darker as compared to ambient teeth. Man was really afraid of preparation. Common treatment with endodontic over post and ceramic crown was refused by patient. Minimal preparation on the porcelain veneer was the only solution.



Dental technic use IPS e.max Press technique. Method of fixing was the combination of the adhesive and the self-adhesive technique.

Porcelain veneer is one of the possibilities how to reach an excellent aesthetics in the frontal part of teeth.

Restoration of the beautiful smile by the porcelain veneers returned both patients self-esteem.